

Matchup Guide

Natrelle® 133 Series Tissue Expanders with Natrelle® 410 Anatomical Gel and Natrelle INSPIRA® Round Cohesive and SoftTouch Gel Breast Implants



ART. SCIENCE. *Natrelle*

First Stage: Natrelle® 133 Series Tissue Expanders With Suture Tabs

Select an expander below based on the base width measurement and find the complementary range of Natrelle® 410 Anatomical Gel Breast Implants to the right in the same color.*

Style 133FV-T[†]
Full Height, Variable Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	11.5	5.0	300	133FV-11-T
12.0	12.5	5.3	400	133FV-12-T
13.0	13.5	5.7	500	133FV-13-T
14.0	14.5	6.2	600	133FV-14-T
15.0	15.5	6.7	750	133FV-15-T
16.0	16.5	6.8	850	133FV-16-T

Style 133FX-T[†]
Full Height, Extra Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	11.5	5.9	350	133FX-11-T
12.0	12.5	6.3	450	133FX-12-T
13.0	13.5	6.7	550	133FX-13-T
14.0	14.5	7.1	650	133FX-14-T
15.0	15.5	7.6	800	133FX-15-T
15.5	16.0	8.1	950	133FX-16-T

Style 133MV-T[†]
Moderate Height, Variable Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	10.0	4.9	250	133MV-11-T
12.0	11.0	5.2	300	133MV-12-T
13.0	12.0	5.6	400	133MV-13-T
14.0	13.0	6.0	500	133MV-14-T
15.0	14.0	6.3	600	133MV-15-T
16.0	15.0	6.6	700	133MV-16-T

Style 133MX-T[†]
Moderate Height, Extra Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	10.0	5.9	300	133MX-11-T
12.0	11.0	6.3	400	133MX-12-T
13.0	12.0	6.7	500	133MX-13-T
14.0	13.0	7.1	600	133MX-14-T
15.0	14.0	7.6	700	133MX-15-T
15.5	14.5	8.1	850	133MX-16-T

Style 133SV-T[†]
Short Height, Variable Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	7.6	4.7	200	133SV-11-T
12.0	8.4	5.3	250	133SV-12-T
13.0	9.1	5.7	300	133SV-13-T
14.0	9.7	6.0	375	133SV-14-T
15.0	10.5	6.3	450	133SV-15-T
16.0	11.3	6.3	550	133SV-16-T

Style 133SX-T[†]
Short Height, Extra Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	9.0	5.9	250	133SX-11-T
12.0	10.0	6.3	350	133SX-12-T
13.0	11.0	6.7	400	133SX-13-T
14.0	12.0	7.1	500	133SX-14-T
15.0	13.0	7.6	650	133SX-15-T
15.5	13.5	8.1	800	133SX-16-T

Style 133LV-T[†]
Low Height, Variable Projection

WIDTH (CM)	HEIGHT (CM) INTRA-ARC	PROJECTION (CM)	HEIGHT (CM) TIP TO BASE	FILL VOLUME (CC)	CATALOG NUMBER
11.0	7.0	4.7	7.6	150	133LV-11-T
12.0	7.5	5.3	8.4	200	133LV-12-T
13.0	8.0	5.7	9.1	300	133LV-13-T
14.0	8.5	6.0	9.7	350	133LV-14-T
15.0	9.0	6.3	10.5	400	133LV-15-T
16.0	9.5	6.3	11.3	500	133LV-16-T



Natrelle® Breast Implants Important Information

INDICATIONS

Natrelle® Breast Implants are indicated for women for the following:

- Breast augmentation for women at least 22 years old for silicone-filled implants.
- Breast augmentation for women at least 18 years old for saline-filled implants.

Breast augmentation includes primary breast augmentation to increase breast size, as well as revision surgery to correct or improve the result of a primary breast augmentation surgery.

- **Breast reconstruction.** Breast reconstruction includes primary reconstruction to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the result of a primary breast reconstruction surgery.

Please see *Natrelle*® 133 Tissue Expanders Indications and Important Safety Information on other side.

Second Stage: Natrelle® 410 Anatomical Gel Breast Implants

These color-coded charts are intended as guides to assist you in selecting complementary expanders and implants.

410 Style FL
Full Height, Low Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.0	10.5	2.9	140	FL-410140
11.0	11.5	3.1	190	FL-410190
11.5	12.0	3.2	220	FL-410220
12.0	12.5	3.4	250	FL-410250
13.0	13.5	3.8	320	FL-410320

410 Style ML
Moderate Height, Low Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.0	9.1	2.9	125	ML-410125
11.0	10.1	3.1	170	ML-410170
11.5	10.6	3.2	195	ML-410195
12.0	11.1	3.4	220	ML-410220
13.0	12.1	3.8	285	ML-410285

410 Style LL
Low Height, Low Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.5	8.6	3.0	135	LL-410135
11.5	9.6	3.2	180	LL-410180
12.0	10.1	3.4	210	LL-410210
12.5	10.5	3.6	240	LL-410240
13.5	11.4	4.0	300	LL-410300

IMPORTANT SAFETY INFORMATION

- CONTRAINDICATIONS**
Breast implant surgery should not be performed in:
- Women with active infection anywhere in their body.
 - Women with existing cancer or precancer of their breast who have not received adequate treatment for those conditions.
 - Women who are currently pregnant or nursing.

- WARNINGS**
- Breast implants are not lifetime devices or necessarily a one-time surgery.
 - Avoid damage during surgery: Care should be taken to avoid the use of excessive force and to minimize handling of the implant. Use care when using surgical instruments in proximity with the breast implant. For more information, please see the full Directions for Use.

- PRECAUTIONS**
Safety and effectiveness have not been established in patients with the following:
- Autoimmune diseases (eg, lupus and scleroderma).
 - A compromised immune system (eg, currently receiving immunosuppressive therapy).
 - Planned chemotherapy following breast implant placement.
 - Planned radiation therapy to the breast following breast implant placement.
 - Conditions or medications that interfere with wound healing and blood clotting.
 - Reduced blood supply to breast tissue.

410 Style FM
Full Height, Moderate Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.5	11.0	3.8	205	FM-410205
11.0	11.5	4.0	235	FM-410235
11.5	12.0	4.2	270	FM-410270
12.0	12.5	4.4	310	FM-410310
12.5	13.0	4.6	350	FM-410350
13.0	13.5	4.8	395	FM-410395
13.5	14.0	5.0	440	FM-410440
14.0	14.5	5.2	500	FM-410500
14.5	15.0	5.4	550	FM-410550
15.0	15.5	5.5	605	FM-410605
15.5	16.0	5.6	670	FM-410670

410 Style MM
Moderate Height, Moderate Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.0	9.1	3.6	160	MM-410160
10.5	9.6	3.8	185	MM-410185
11.0	10.1	4.0	215	MM-410215
11.5	10.6	4.2	245	MM-410245
12.0	11.1	4.4	280	MM-410280
12.5	11.6	4.6	320	MM-410320
13.0	12.1	4.8	360	MM-410360
13.5	12.5	5.0	400	MM-410400
14.0	12.9	5.2	450	MM-410450

410 Style LM
Low Height, Moderate Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.0	8.1	3.6	140	LM-410140
11.0	9.1	4.0	190	LM-410190
11.5	9.6	4.2	220	LM-410220
12.0	10.1	4.4	250	LM-410250
13.0	10.9	4.8	320	LM-410320

410 Style MF
Moderate Height, Full Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
9.5	8.6	3.7	140	MF-410140
10.0	9.1	4.0	165	MF-410165
10.5	9.6	4.2	195	MF-410195
11.0	10.1	4.4	225	MF-410225
11.5	10.6	4.6	255	MF-410255
12.0	11.1	4.8	295	MF-410295
12.5	11.6	5.1	335	MF-410335
13.0	12.1	5.2	375	MF-410375
13.5	12.5	5.3	420	MF-410420
14.0	12.9	5.6	470	MF-410470
14.5	13.2	5.8	525	MF-410525
15.0	13.6	6.1	580	MF-410580
15.5	13.9	6.2	640	MF-410640

410 Style LF
Low Height, Full Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
9.5	7.6	3.7	125	LF-410125
10.0	8.1	4.0	150	LF-410150
10.5	8.6	4.2	175	LF-410175
11.0	9.1	4.4	205	LF-410205
11.5	9.6	4.6	240	LF-410240
12.0	10.1	4.8	270	LF-410270
12.5	10.5	5.1	310	LF-410310
13.5	11.4	5.3	390	LF-410390
14.0	11.8	5.6	440	LF-410440
14.5	12.2	5.8	490	LF-410490
15.0	12.6	6.1	540	LF-410540
15.5	13.0	6.2	595	LF-410595

Other potential adverse events that may occur with breast implant surgery include: asymmetry, breast pain, breast/skin sensation changes, capsular calcification, delayed wound healing, hematoma, hypertrophic scarring/scarring, implant extrusion, implant malposition, implant palpability/visibility, infection, nipple complications, redness, seroma, swelling, tissue/skin necrosis, wrinkling/rippling.

410 Style FF
Full Height, Full Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.0	10.5	4.0	185	FF-410185
10.5	11.0	4.2	220	FF-410220
11.0	11.5	4.4	255	FF-410255
11.5	12.0	4.6	290	FF-410290
12.0	12.5	4.8	335	FF-410335
12.5	13.0	5.1	375	FF-410375
13.0	13.5	5.2	425	FF-410425
13.5	14.0	5.3	475	FF-410475
14.0	14.5	5.6	535	FF-410535
14.5	15.0	5.8	595	FF-410595
15.0	15.5	6.1	655	FF-410655
15.5	16.0	6.2	740	FF-410740

410 Style MM
Moderate Height, Moderate Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
10.0	9.1	3.6	160	MM-410160
10.5	9.6	3.8	185	MM-410185
11.0	10.1	4.0	215	MM-410215
11.5	10.6	4.2	245	MM-410245
12.0	11.1	4.4	280	MM-410280
12.5	11.6	4.6	320	MM-410320
13.0	12.1	4.8	360	MM-410360
13.5	12.5	5.0	400	MM-410400
14.0	12.9	5.2	450	MM-410450

410 Style MF
Moderate Height, Full Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
9.5	8.6	3.7	140	MF-410140
10.0	9.1	4.0	165	MF-410165
10.5	9.6	4.2	195	MF-410195
11.0	10.1	4.4	225	MF-410225
11.5	10.6	4.6	255	MF-410255
12.0	11.1	4.8	295	MF-410295
12.5	11.6	5.1	335	MF-410335
13.0	12.1	5.2	375	MF-410375
13.5	12.5	5.3	420	MF-410420
14.0	12.9	5.6	470	MF-410470
14.5	13.2	5.8	525	MF-410525
15.0	13.6	6.1	580	MF-410580
15.5	13.9	6.2	640	MF-410640

410 Style LF
Low Height, Full Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
9.5	7.6	3.7	125	LF-410125
10.0	8.1	4.0	150	LF-410150
10.5	8.6	4.2	175	LF-410175
11.0	9.1	4.4	205	LF-410205
11.5	9.6	4.6	240	LF-410240
12.0	10.1	4.8	270	LF-410270
12.5	10.5	5.1	310	LF-410310
13.5	11.4	5.3	390	LF-410390
14.0	11.8	5.6	440	LF-410440
14.5	12.2	5.8	490	LF-410490
15.0	12.6	6.1	540	LF-410540
15.5	13.0	6.2	595	LF-410595


410 Style FX
Full Height, Extra-Full Projection

WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	CATALOG NUMBER
9.5	10.0	4.6	185	FX-410185
10.0	10.5	4.9	215	FX-4102

First Stage: Natrelle® 133 Series Tissue Expanders With Suture Tabs


Select an expander below based on the base width measurement and find the complementary range of Natrelle INSPIRA® Round Gel Breast Implants to the right in the same color.

Style 133FV-T*
Full Height, Variable Projection




WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	11.5	5.0	300	133FV-11-T
12.0	12.5	5.3	400	133FV-12-T
13.0	13.5	5.7	500	133FV-13-T
14.0	14.5	6.2	600	133FV-14-T
15.0	15.5	6.7	750	133FV-15-T
16.0	16.5	6.8	850	133FV-16-T

Style 133MV-T*
Moderate Height, Variable Projection




WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	10.0	4.9	250	133MV-11-T
12.0	11.0	5.2	300	133MV-12-T
13.0	12.0	5.6	400	133MV-13-T
14.0	13.0	6.0	500	133MV-14-T
15.0	14.0	6.3	600	133MV-15-T
16.0	15.0	6.6	700	133MV-16-T

Style 133SV-T*
Short Height, Variable Projection




WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	7.6	4.7	200	133SV-11-T
12.0	8.4	5.3	250	133SV-12-T
13.0	9.1	5.7	300	133SV-13-T
14.0	9.7	6.0	375	133SV-14-T
15.0	10.5	6.3	450	133SV-15-T
16.0	11.3	6.3	550	133SV-16-T

Style 133LV-T*
Low Height, Variable Projection




WIDTH (CM)	HEIGHT (CM) INTRA-ARC	PROJECTION (CM)	HEIGHT (CM) TIP TO BASE	FILL VOLUME (CC)	CATALOG NUMBER
11.0	7.0	4.7	7.6	150	133LV-11-T
12.0	7.5	5.3	8.4	200	133LV-12-T
13.0	8.0	5.7	9.1	300	133LV-13-T
14.0	8.5	6.0	9.7	350	133LV-14-T
15.0	9.0	6.3	10.5	400	133LV-15-T
16.0	9.5	6.3	11.3	500	133LV-16-T

Style 133FX-T*
Full Height, Extra Projection




WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	11.5	5.9	350	133FX-11-T
12.0	12.5	6.3	450	133FX-12-T
13.0	13.5	6.7	550	133FX-13-T
14.0	14.5	7.1	650	133FX-14-T
15.0	15.5	7.6	800	133FX-15-T
15.5	16.0	8.1	950	133FX-16-T

Style 133MX-T*
Moderate Height, Extra Projection



WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	10.0	5.9	300	133MX-11-T
12.0	11.0	6.3	400	133MX-12-T
13.0	12.0	6.7	500	133MX-13-T
14.0	13.0	7.1	600	133MX-14-T
15.0	14.0	7.6	700	133MX-15-T
15.5	14.5	8.1	850	133MX-16-T

Style 133SX-T*
Short Height, Extra Projection




WIDTH (CM)	HEIGHT (CM)	PROJECTION (CM)	FILL VOLUME (CC)	CATALOG NUMBER
11.0	9.0	5.9	250	133SX-11-T
12.0	10.0	6.3	350	133SX-12-T
13.0	11.0	6.7	400	133SX-13-T
14.0	12.0	7.1	500	133SX-14-T
15.0	13.0	7.6	650	133SX-15-T
15.5	13.5	8.1	800	133SX-16-T



Second Stage: Natrelle INSPIRA® Smooth Round Cohesive and SoftTouch Gel Breast Implants


These color-coded charts are intended as guides to assist you in selecting complementary expanders and implants.

Natrelle INSPIRA® Styles SCM & SSM
Smooth, Moderate Profile




DIAMETER (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	Smooth Cohesive (SCM)	Smooth SoftTouch (SSM)
9.50	3.0	140	SCM-140	SSM-140
10.00	3.3	175	SCM-175	SSM-175
10.25	3.4	195	SCM-195	SSM-195
10.50	3.5	210	SCM-210	SSM-210
11.00	3.7	240	SCM-240	SSM-240
11.25	3.8	255	SCM-255	SSM-255
11.50	3.9	275	SCM-275	SSM-275
11.75	3.9	295	SCM-295	SSM-295
12.00	4.0	310	SCM-310	SSM-310
12.25	4.1	330	SCM-330	SSM-330
12.50	4.2	345	SCM-345	SSM-345
12.75	4.2	360	SCM-360	SSM-360
13.00	4.3	375	SCM-375	SSM-375
13.25	4.4	405	SCM-405	SSM-405
13.50	4.5	445	SCM-445	SSM-445
14.00	4.6	485	SCM-485	SSM-485
14.50	4.8	520	SCM-520	SSM-520
14.75	4.9	560	SCM-560	SSM-560
15.00	4.9	600	SCM-600	SSM-600
15.50	5.0	640	SCM-640	SSM-640
16.00	5.2	685	SCM-685	SSM-685
16.50	5.0	755	SCM-755	SSM-755

Natrelle INSPIRA® Styles SCL & SSL
Smooth, Low Profile




DIAMETER (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	Smooth Cohesive (SCL)	Smooth SoftTouch (SSL)
10.00	2.0	110	SCL-110	SSL-110
10.25	2.0	125	SCL-125	SSL-125
10.50	2.1	140	SCL-140	SSL-140
11.00	2.2	170	SCL-170	SSL-170
11.50	2.3	200	SCL-200	SSL-200
12.00	2.4	230	SCL-230	SSL-230
12.50	2.5	260	SCL-260	SSL-260
13.00	2.6	290	SCL-290	SSL-290
13.50	2.7	320	SCL-320	SSL-320
14.00	2.8	350	SCL-350	SSL-350
14.50	2.9	380	SCL-380	SSL-380
15.00	3.0	410	SCL-410	SSL-410
15.50	3.1	460	SCL-460	SSL-460
16.00	3.2	510	SCL-510	SSL-510
16.50	3.3	560	SCL-560	SSL-560
17.00	3.4	610	SCL-610	SSL-610

Natrelle INSPIRA® Styles SCF & SSF
Smooth, Full Profile




DIAMETER (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	Smooth Cohesive (SCF)	Smooth SoftTouch (SSF)
9.50	4.0	180	SCF-180	SSF-180
9.75	4.1	200	SCF-200	SSF-200
10.00	4.2	220	SCF-220	SSF-220
10.50	4.3	240	SCF-240	SSF-240
11.00	4.4	265	SCF-265	SSF-265
11.25	4.5	295	SCF-295	SSF-295
11.40	4.6	325	SCF-325	SSF-325
11.50	4.7	335	SCF-335	SSF-335
11.75	4.8	345	SCF-345	SSF-345
12.00	4.9	365	SCF-365	SSF-365
12.25	5.0	385	SCF-385	SSF-385
12.50	5.1	415	SCF-415	SSF-415
12.75	5.3	450	SCF-450	SSF-450
13.00	5.4	485	SCF-485	SSF-485
13.25	5.5	520	SCF-520	SSF-520
13.50	5.7	560	SCF-560	SSF-560
14.00	5.8	605	SCF-605	SSF-605
14.50	5.9	650	SCF-650	SSF-650
14.75	6.0	695	SCF-695	SSF-695
15.00	6.1	745	SCF-745	SSF-745
15.50	6.0	770	SCF-770	SSF-770

Natrelle INSPIRA® Styles SCLP & SSLP
Smooth, Low-Plus Profile



DIAMETER (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	Smooth Cohesive (SCLP)	Smooth SoftTouch (SSLP)
9.50	2.5	125	SCLP-125	SSLP-125
10.00	2.6	145	SCLP-145	SSLP-145
10.50	2.7	165	SCLP-165	SSLP-165
11.00	2.8	190	SCLP-190	SSLP-190
11.25	2.9	205	SCLP-205	SSLP-205
11.50	3.0	220	SCLP-220	SSLP-220
11.75	3.1	235	SCLP-235	SSLP-235
12.00	3.1	250	SCLP-250	SSLP-250
12.25	3.2	265	SCLP-265	SSLP-265
12.50	3.3	280	SCLP-280	SSLP-280
12.75	3.3	300	SCLP-300	SSLP-300
13.00	3.4	320	SCLP-320	SSLP-320
13.25	3.5	340	SCLP-340	SSLP-340
13.50	3.5	360	SCLP-360	SSLP-360
14.00	3.6	400	SCLP-400	SSLP-400
14.50	3.8	440	SCLP-440	SSLP-440
15.00	3.9	490	SCLP-490	SSLP-490
15.50	4.0	540	SCLP-540	SSLP-540
16.00	4.2	590	SCLP-590	SSLP-590
16.50	4.3	640	SCLP-640	SSLP-640

Natrelle INSPIRA® Styles SCX & SSX
Smooth, Extra-Full Profile



DIAMETER (CM)	PROJECTION (CM)	IMPLANT VOLUME (CC)	Smooth Cohesive (SCX)	Smooth SoftTouch (SSX)
9.50	4.2	200	SCX-200	SSX-200
10.00	4.5	225	SCX-225	SSX-225
10.25	4.8	255	SCX-255	SSX-255
10.50	5.0	285	SCX-285	SSX-285
10.75	5.2	310	SCX-310	SSX-310
11.00	5.4	340	SCX-340	SSX-340
11.25	5.6	375	SCX-375	SSX-375
11.50	5.8	400	SCX-400	SSX-400
11.75	5.9	420	SCX-420	SSX-420
12.00	6.0	445	SCX-445	SSX-445
12.25	6.1	470	SCX-470	SSX-470
12.50	6.1	495	SCX-495	SSX-495
12.75	6.2	525	SCX-525	SSX-525
13.00	6.2	545	SCX-545	SSX-545
13.25	6.3	560	SCX-560	SSX-560
13.50	6.4	580	SCX-580	SSX-580
13.75	6.4	615	SCX-615	SSX-615
14.00	6.5	650	SCX-650	SSX-650
14.25	6.6	700	SCX-700	SSX-700
14.50	6.7	750	SCX-750	SSX-750
15.50	6.7	800	SCX-800	SSX-800

*Each style of Natrelle® 133 Tissue Expanders is available with or without suture tabs.

Also available in textured styles!

Please see the Natrelle® Product Catalog to view Natrelle INSPIRA® Textured Round Gel Breast Implants.

IMPORTANT SAFETY INFORMATION (continued) WARNINGS (continued)

- DO NOT alter the tissue expander or use adulterated fill.
- Fill only with sterile saline for injection as described in INSTRUCTIONS FOR USE. DO NOT expose to contaminants.
- DO NOT expand if the pressure will compromise wound healing or vasculature of overlying tissue, or beyond patient or tissue tolerance. Stop filling immediately if tissue damage, wound dehiscence, abnormal skin pallor, erythema, edema, pain, or tenderness are observed.
- Active infection anywhere may increase risk of periprosthetic infection. Postoperative infections should be treated aggressively. Unresponsive or necrotizing infection may require premature removal.
- Natrelle® 133 Tissue Expanders are temporary, single-use only devices, and are not to be used for permanent implantation or beyond 6 months. Tissue expansion in breast reconstruction typically requires 4 months to 6 months.
- When using suturing tabs be careful to avoid piercing the shell. Use a new one if damage occurs.

PRECAUTIONS

Active infections may need to be treated and resolved before surgery. Allergan relies on the surgeon to know and follow proper surgical procedures and carefully evaluate patient suitability using standard practice and individual experience. Avoid damage to the tissue expander and use a sterile backup in case of damage. Pay careful attention to tissue tolerance and hemostasis during surgery. Expansion should proceed moderately and never beyond patient or tissue tolerance. Avoid contamination in any postoperative procedure.

ADVERSE REACTIONS

Deflation, tissue damage, infection, extrusion, hematoma/seroma, capsular contracture, premature explantation, displacement, effects on bone, pain, sensation, distortion, inadequate tissue flap, and inflammatory reaction.

For more information, please visit www.allergan.com/labeling/usa.htm. To report a problem with Natrelle®, please call Allergan at 1-800-433-8871.

Natrelle® 133 Tissue Expanders are available by prescription only.

WARNINGS

- DO NOT use Natrelle® 133 Tissue Expanders in patients who already have implanted devices that would be affected by a magnetic field (see Contraindications), because the MAGNA-SITE® integrated injection site contains a strong rare-earth, permanent magnet. Diagnostic testing with Magnetic Resonance Imaging (MRI) is contraindicated in patients with Natrelle® 133 Tissue Expanders in place.



Natrelle® 133 Tissue Expanders With/Without Suture Tabs and With MAGNA-SITE® Injection Sites Important Information

INDICATIONS

Natrelle® 133 Tissue Expanders are indicated for:

- Breast reconstruction following mastectomy.
- Treatment of underdeveloped breasts.
- Treatment of soft tissue deformities.

IMPORTANT SAFETY INFORMATION CONTRAINDICATIONS

Natrelle® 133 Tissue Expanders SHOULD NOT be used in patients:

- Who already have implanted devices that would be affected by a magnetic field (eg, pacemakers, drug infusion devices, artificial sensing devices).

- Whose tissue at the expansion site is determined to be unsuitable.
- Who have an active infection or a residual gross tumor at the expansion site.
- Undergoing adjuvant radiation therapy.

- Whose physiological condition (eg, sensitive over- or underlying anatomy, obesity, smoking, diabetes, autoimmune disease, hypertension, chronic lung or severe cardiovascular disease, or osteogenesis imperfecta) or use of certain drugs (including those that interfere with blood clotting or affect tissue viability) poses an unduly high risk of surgical and/or postoperative complications.
- Who are psychologically unsuitable.

Please see Natrelle® Breast Implants Indications and Important Safety Information on other side.

